



ZETIN

# Glazing Furnace

ZTCF-20B Mini



## Technical Parameters

Max.Temperature	1100°C
Temperature Accuracy	Within $\pm 1^{\circ}\text{C}$
Power	700W
Voltage	220V 50/60HZ
Chamber Size	60×70×60mm(L×W×H)
Size	220×170×266mm(L×W×H)
Net Weight	8.5 KG



## Features

### High-Quality and Efficiency

Imported heating wire, uniform heating, long service life, and pollution-free.

### Precision Temperature Control

PID intelligent temperature control system with high temperature control accuracy, of  $\pm 1^{\circ}\text{C}$ .

### Program Preset Function

100 programs can be stored and retrieved easily.

### High-Definition Touch Screen

Large touch screen, parameterized design interface, and panoramic display of all parameters.

### Energy-Saving

New environmentally friendly insulation materials with good insulation performance ensure stable furnace temperature.

### Small Size

Compact and lightweight, suitable for chairside.